



PBL InterferonSource

MOUSE MONOCLONAL ANTIBODY AGAINST HUMAN INTERFERON GAMMA RECEPTOR (MMHG-2)

Product Number: 21550-1

Lot Number:

Size: 100 µg

Product Information

Description: Mouse Anti-Human Interferon Gamma, non-neutralizing

Clone: MMHG-2

Contents: 100 µg / vial

Volume: 0.1 ml

Concentration: 1 mg/ml

Buffer: phosphate buffered saline (PBS) containing 0.1% bovine serum albumin (BSA)

Antigen: human interferon gamma

Isotype: IgG₁

Bioactivity: binds to but does not neutralize human interferon gamma

Assay Used to Measure Bioactivity: One neutralization unit is the amount of antiserum required to neutralize one unit of human interferon gamma (Hu-IFN- γ) to a 50% endpoint. Interferon was titrated with the use of the cytopathic effect inhibition assay as described [Rubinstein, S., Familletti, P.C., and Pestka, S. (1981) "Convenient Assay for Interferons," *J. Virol.* **37**, 755-758; Familletti, P.C., Rubinstein, S., and Pestka, S. (1981) "A Convenient and Rapid Cytopathic Effect Inhibition Assay for Interferon," in *Methods in Enzymology*, Vol. 78 (S. Pestka, ed.), Academic Press, New York, 387-394]. In this antiviral assay for interferon about 1 unit/ml of interferon is the quantity necessary to produce a cytopathic effect of 50%. The units are determined with respect to the international reference standard for Hu-IFN- γ provided by the National Institutes of Health [see Pestka, S. (1986) "Interferon Standards and General Abbreviations," in *Methods in Enzymology* (S. Pestka, ed.), Academic Press, New York **119**, 14-23].

Suggested Applications: Western blot; immunoprecipitation; immunohistochemistry; ELISA

Please note that these applications are presented for suggested use only and have not been fully evaluated by PBL

Shipping Information

Shipping Conditions: Wet ice

Physical State of Product During Shipping: Liquid

Special Conditions/Comments: After receipt, this product should be kept at -20°C for retention of full activity. When thawing, the contents of the tube should be apportioned in separate tubes so that freezing and thawing is kept to a minimum.

**This Product is for RESEARCH USE ONLY and is not for sale or use in any commercial kit.
NOT FOR DIAGNOSTIC OR THERAPEUTIC USE, OR USE IN HUMANS.**

v.1.2

your total interferon source

131 Ethel Road West, Suite 6 • Piscataway, NJ 08854 USA • Telephone: 732.777.9123 • Facsimile: 732.777.9141
Toll Free: 877.PBL.8881 • email: info@interferonsource.com